



# Steering devices

Driving without limitations

BEVER MOBILITY PRODUCTS



# Steering devices

## Driving without limitations

The various types of steering devices allow safe operation of the steering wheel for one-handed drivers, so they can remain in control of the steering at all times. Abled bodied drivers, can easily remove the steering device with the quick release button.



For the driver with full use of hand and arm when steering. The steering knob device is designed to fit comfortably in the palm of the hand, giving the driver extra power and control when operating the steering wheel.



The 1-pin device is an alternative to the steering knob, for the driver who cannot fully open his hand. The padded disc provides comfortable hand support.

# Steering devices

## Driving without limitations



The V-grip device is for the driver with minimal grip and reduced hand function. It can be adjusted to the size of the hand. The hand remains securely and comfortably in place.



The 3-pin device is for the driver with adequate arm strength, but minimal grip strength and / or reduced wrist stability. It supports the hand and wrist and keeps them comfortably in position. The tri-pin device consists of 3 legs. The front leg has a fixed position, the rear legs can be adjusted inward, outward, forward and backward. This allows the tri-pin steering grip to be perfectly adjustable to fit the shape of the hand.



The palm grip device is an alternative to the 3-pin device, if the driver has no grip and / or reduced wrist stability. The hand is stabilized by a cuff.

# Steering devices

## Driving without limitations

- Ergonomic design
- Easy to remove with a quick release button
- Fits most steering wheel sizes
- Fits multiple positions on steering wheel
- Same steering-adapter as with SmartSteer



Complete your steering devices with SmartSteer



**Bever Mobility Products**  
2885 Sanford Ave, SW# 41693  
Grandville, MI 49418  
**Tel. (888) 959 6198**